

Product datasheet

info@arigobio.com

ARG80437 Human 3 CAT (Adrenaline + Noradrenaline + Dopamine) ELISA Kit

Package: 96 wells Store at: 4°C

Summary

Product Description ARG80437 Human 3 CAT (Adrenaline + Noradrenaline + Dopamine) ELISA Kit is an Enzyme

Immunoassay kit for the quantification of Adrenaline (Epinephrine), Noradrenaline (Norepinephrine),

and Dopamine in plasma (EDTA) and urine.

This kit includes 96 wells for Adrenaline, 96 wells for Noradrenaline and 96 wells for Dopamine.

Tested Reactivity Hu

Tested Application ELISA

Target Name 3 CAT (Adrenaline + Noradrenaline + Dopamine)

Conjugation HRP

Sensitivity AD: 0.8 ng/ml (urine); 0.01 pg/ml (plasma).

NAD: 1.7 ng/ml (urine); 0.04 pg/ml (plasma). DOP: 2.5 ng/ml (urine); 1.9 pg/ml (plasma).

Sample Type Plasma (EDTA) and urine.

Standard Range AD: 1 - 200 ng/ml

NAD: 5 - 1000 ng/ml DOP: 4.5 - 2000 ng/ml

Sample Volume 10 µl urine; 300 µl plasma

Application Instructions

Assay Time 30 min, 2 h, 30, 25 min (RT/shaker)

Properties

Form 96 well

Storage instruction Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test

reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual

for detail temperatures of the components.

Note For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight Related products:


3 CAT ELISA Kits; Related news:

Astrocyte-to-neuron conversion for Parkinson's disease treatment

New ELISA data calculation tool: Simplify the ELISA analysis by GainData


Research Area Neuroscience kit; Signaling Transduction kit

www.arigobio.com arigo.nuts about antibodies 1/2




ARG80437 Human 3 CAT (Adrenaline + Noradrenaline + Dopamine) ELISA Kit standard curve example image for Adrenaline

The representative standard curve of ARG80437 Human 3 CAT (Adrenaline + Noradrenaline + Dopamine) ELISA Kit. The standard curve is for demonstration only and cannot be used in place of data generations at the time of assay. The standard curve should be generated each time the assay is performed.



ARG80437 Human 3 CAT (Adrenaline + Noradrenaline + Dopamine) ELISA Kit standard curve example image for Noradrenaline

The representative standard curve of ARG80437 Human 3 CAT (Adrenaline + Noradrenaline + Dopamine) ELISA Kit. The standard curve is for demonstration only and cannot be used in place of data generations at the time of assay. The standard curve should be generated each time the assay is performed.



ARG80437 Human 3 CAT (Adrenaline + Noradrenaline + Dopamine) ELISA Kit standard curve example image for Dopamine

The representative standard curve of ARG80437 Human 3 CAT (Adrenaline + Noradrenaline + Dopamine) ELISA Kit. The standard curve is for demonstration only and cannot be used in place of data generations at the time of assay. The standard curve should be generated each time the assay is performed.