

Product datasheet

info@arigobio.com

ARG83003 Arg8 Vasopressin ELISA Kit

Package: 96 wells Store at: 4°C

Summary

Product Description ARG83003 Arg8 Vasopressin ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human,

Mouse, Rat, Sheep, Dog Arg8 Vasopressin in plasma.

Tested Reactivity Hu, Ms, Rat, Dog, Sheep

Tested Application ELISA

Specificity Cross-Reactivity:

[Arg8]-Vasopressin (Human, Mouse, Rat, Ovine, Canine): 100%; [Arg8]-Vasotocin: 54%;

[Lys8]-Vasopressin (Human, Rat, Mouse, Ovine), Deamino-[D-Arg8]-Vassopressin, Oxytocin, LH-RH, ACTH (Human), alpha-ANP (1-28) (Human, Canine), Met-Enkephalin, PACAP38 (Human, Ovine, Rat) and

PACAP27-NH2 (Human, Ovine, Rat): 0%.

Target Name Arg8 Vasopressin

Conjugation HRP

Conjugation Note Substrate: TMB and read at 450 nm.

Sensitivity 0.04 ng/ml

Detection Range 0.04 - 0.51 ng/ml

Sample Type Plasma (Sample extraction is recommended)

Standard Range 0.01 - 100 ng/ml

Sample Volume $50 \, \mu l$

Precision Intra-assay: less than 10%

Inter-assay: < 15%

Alternate Names Arginine-vasopressin; Vasopressin-neurophysin 2-Arg8 Vasopressin; Neurophysin-II; AVP; Arginine

Vasopressin; ARVP; ADH; Prepro-Arginine-Vasopressin-Neurophysin II; Vasopressin-Neurophysin 2-Copeptin; Antidiuretic Hormone; Neurohypophyseal; Prepro-AVP-NP; Copeptin; AVP-NPII; VP;

Vasopressin-Neurophysin II-Copeptin; Diabetes Insipidus; Neurophysin II; AVRP

Application Instructions

Assay Time ~ 4 hours

Properties

Form 96 well

Storage instruction Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test

reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual

for detail temperatures of the components.

Note For laboratory research only, not for drug, diagnostic or other use. This kit cannot ship or deliver to

addresses in the following countries: Argentina, Austria, Belgium, China, France, Germany, Hungary, India, Italy, Japan, Liechtenstein, Luxembourg, Mexico, Spain, Switzerland, United Kingdom. Due to the assay method performed by this product is covered by Patent No. EP1738178 in these countries.

www.arigobio.com arigo.nuts about antibodies 1/2

Bioinformation

Gene Symbol

AVP

Gene Full Name

arginine vasopressin

Background

This gene encodes a precursor protein consisting of arginine vasopressin and two associated proteins, neurophysin 2 and a glycopeptide, Arg8 Vasopressin. Arginine vasopressin is a posterior pituitary hormone which is synthesized in the supraoptic nucleus and paraventricular nucleus of the hypothalamus. Along with its carrier protein, neurophysin 2, it is packaged into neurosecretory vesicles and transported axonally to the nerve endings in the neurohypophysis where it is either stored or secreted into the bloodstream. The precursor is thought to be activated while it is being transported along the axon to the posterior pituitary. Arginine vasopressin acts as a growth factor by enhancing pH regulation through acid-base transport systems. It has a direct antidiuretic action on the kidney, and also causes vasoconstriction of the peripheral vessels. This hormone can contract smooth muscle during parturition and lactation. It is also involved in cognition, tolerance, adaptation and complex sexual and maternal behaviour, as well as in the regulation of water excretion and cardiovascular functions. Mutations in this gene cause autosomal dominant neurohypophyseal diabetes insipidus (ADNDI). [provided by RefSeq, Mar 2010]

Function

Arg-vasopressin has a direct antidiuretic action on the kidney, it also causes vasoconstriction of the peripheral vessels. Acts by binding to vasopressin receptors (V1bR/AVPR1B, V1aR/AVPR1A, and V2R/AVPR2) (PubMed:18174156) [UniProt]


Highlight

Related products:

<u>Vasopressin antibodies;</u> <u>Vasopressin ELISA Kits;</u> <u>Vasopressin Duos / Panels;</u>

New ELISA data calculation tool: Simplify the ELISA analysis by GainData

Images



ARG83003 Arg8 Vasopressin ELISA Kit standard curve image

ARG83003 Arg8 Vasopressin ELISA Kit results of a typical standard run with optical density reading at 450 nm.