

ARG10063 anti-HCV Core antigen / HCcAg antibody [A1/3D1]

Package: 100 µg

Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [A1/3D1] recognizes Hepatitis C virus core antigen
Tested Reactivity	HCV
Tested Application	ELISA
Specificity	No cross reaction with recombinant or synthetic HCV non-structural proteins (NS-3 and NS-4).
Host	Mouse
Clonality	Monoclonal
Clone	A1/3D1
Isotype	IgG1, kappa
Target Name	HCV Core antigen / HCcAg
Species	HCV
Immunogen	Synthetic peptides derived from HCV core protein
Conjugation	Un-conjugated

Application Instructions

Application Note