

ARG10142
anti-FK506 antibody [FK1]Package: 100 µg, 50 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [FK1] recognizes FK506
Tested Reactivity	Other
Tested Application	ELISA
Host	Mouse
Clonality	Monoclonal
Clone	FK1
Isotype	IgM
Target Name	FK506
Species	Others
Immunogen	BSA-conjugated FK-506
Conjugation	Un-conjugated

Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Euglobulin precipitation and dialysis.
Buffer	50 mM Tris-HCl (pH 8.0), 150 mM NaCl and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	1.0-2.0 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Immune System antibody
---------------	------------------------