

ARG10229 anti-Cortisol antibody [CORT-1]

Package: 500 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [CORT-1] recognizes Cortisol
Tested Reactivity	Hu
Tested Application	ELISA
Host	Mouse
Clonality	Monoclonal
Clone	CORT-1
Isotype	IgG1
Target Name	Cortisol
Species	Human
Immunogen	conjugate of purified cortisol with BSA
Conjugation	Un-conjugated

Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Protein G affinity purified
Buffer	PBS (pH 7.4) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	1.0-2.0 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Metabolism antibody; Signaling Transduction antibody
---------------	--