

ARG10267 anti-Staphylococcus Enterotoxin B antibody [S222]

Package: 100 µg, 50 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [S222] recognizes Staphylococcus Enterotoxin B
Tested Reactivity	Bacteria
Tested Application	ELISA
Specificity	Do not cross-react with St. aureus enterotoxins A, C, D and E.
Host	Mouse
Clonality	Monoclonal
Clone	S222
Isotype	IgG1
Target Name	Staphylococcus Enterotoxin B
Immunogen	recombinant St. aureus enterotoxin B expressed in E. coli
Conjugation	Un-conjugated

Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Protein G affinity purified
Buffer	PBS (pH 7.4) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	1.0-2.0 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Microbiology and Infectious Disease antibody
---------------	--