

## ARG10336 anti-Streptavidin antibody [S8E4]

Package: 100 µg, 50 µg  
Store at: -20°C

### Summary

|                     |                                                          |
|---------------------|----------------------------------------------------------|
| Product Description | Mouse Monoclonal antibody [S8E4] recognizes Streptavidin |
| Tested Reactivity   | Bacteria                                                 |
| Tested Application  | ELISA, FACS, IHC-Fr, IHC-P, WB                           |
| Specificity         | No cross-reactivity with egg white avidin                |
| Host                | Mouse                                                    |
| Clonality           | Monoclonal                                               |
| Clone               | S8E4                                                     |
| Isotype             | IgG1                                                     |
| Target Name         | Streptavidin                                             |
| Species             | Bacteria                                                 |
| Immunogen           | Streptavidin                                             |
| Conjugation         | Un-conjugated                                            |

### Application Instructions

|                  |                                                                                                                                            |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Application Note | * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist. |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------|

### Properties

|                     |                                                                                                                                                                                                                                                                                                                            |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form                | Liquid                                                                                                                                                                                                                                                                                                                     |
| Purification        | Protein A affinity purified.                                                                                                                                                                                                                                                                                               |
| Buffer              | PBS (pH 7.4) and 0.1% Sodium azide                                                                                                                                                                                                                                                                                         |
| Preservative        | 0.1% Sodium azide                                                                                                                                                                                                                                                                                                          |
| Concentration       | 1.0-2.0 mg/ml                                                                                                                                                                                                                                                                                                              |
| Storage instruction | For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use. |
| Note                | For laboratory research only, not for drug, diagnostic or other use.                                                                                                                                                                                                                                                       |

### Bioinformation

|                |              |
|----------------|--------------|
| Gene Symbol    | BMA3378      |
| Gene Full Name | streptavidin |

