

ARG10629 anti-Influenza A Haemagglutinin H7 antibody [InA334]

Package: 100 µg, 50 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [InA334] recognizes Influenza A Haemagglutinin H7
Tested Reactivity	Influenza A
Tested Application	ELISA, FLISA, I-ELISA
Specificity	Influenza A haemagglutinin H7
Host	Mouse
Clonality	Monoclonal
Clone	InA334
Isotype	IgG1
Target Name	Influenza A Haemagglutinin H7
Species	Virus
Immunogen	Influenza A haemagglutinin H7.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FLISA	Assay-dependent
	I-ELISA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purification with Protein A.
Buffer	PBS (pH 7.4) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Highlight


Related Antibody Duos and Panels:

[ARG30287 Influenza A Haemagglutinin H7 ELISA Antibody Duo](#)

Related products:


[Influenza A antibodies](#); [Influenza A ELISA Kits](#); [Influenza A Duos / Panels](#); [Anti-Mouse IgG secondary antibodies](#);

Images




ARG10629 anti-Influenza A Haemagglutinin H7 antibody [InA334]
FLISA image

FLISA: Calibration curve for Influenza A haemagglutinin H7 in sandwich fluoroimmunoassay. Capture MAb: ARG10629 anti-Influenza A Haemagglutinin H7 antibody [InA334] - 1 µg/well. Detection MAb: (Eu-chelate labeled): ARG10628 anti-Influenza A Haemagglutinin H7 antibody [InA331] - 0.1 µg/well. Incubation time: 45 min. Antigens: Influenza A/Netherlands/219/03 H7N7 virus.



ARG10629 anti-Influenza A Haemagglutinin H7 antibody [InA334]
FLISA image

FLISA: Calibration curve for Influenza A haemagglutinin H7 in sandwich fluoroimmunoassay. Capture MAb: ARG10629 anti-Influenza A Haemagglutinin H7 antibody [InA334] - 1 µg/well. Detection MAb: (Eu-chelate labeled): ARG10628 anti-Influenza A Haemagglutinin H7 antibody [InA331] - 0.1 µg/well. Incubation time: 45 min. Antigens: Recombinant H7 (A/Chicken/Netherlands/1/03) (aa 17-527).



ARG10629 anti-Influenza A Haemagglutinin H7 antibody [InA334]
indirect ELISA image

Indirect ELISA: Titration curve of ARG10629 anti-Influenza A Haemagglutinin H7 antibody [InA334] specific to haemagglutinin H7. Antigens: Viruses - 200 ng/well and Recombinant H7 (A/Chicken/Netherlands/1/03)(aa 17-527) - 10 ng/well.