

ARG10632 anti-Influenza B Haemagglutinin antibody [InB190]

Package: 100 µg, 50 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [InB190] recognizes Influenza B Haemagglutinin
Tested Reactivity	Influenza B virus
Tested Application	ELISA, I-ELISA, WB
Specificity	Influenza B haemagglutinin 2
Host	Mouse
Clonality	Monoclonal
Clone	InB190
Isotype	IgG2b
Target Name	Influenza B Haemagglutinin
Species	Virus
Immunogen	Influenza B haemagglutinin 2.
Conjugation	Un-conjugated

Application Instructions


Application table	Application	Dilution
	ELISA	Assay-dependent
	I-ELISA	Assay-dependent
	WB	3 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purification with Protein A.
Buffer	PBS (pH 7.4) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.


Gene Symbol	HA
Gene Full Name	Haemagglutinin

Images



ARG10632 anti-Influenza B Haemagglutinin antibody [InB190] WB image

Western blot: Influenza B HA2 immunodetection. 1 µg of 1) Influenza B/Tokio/53/99, 2) Influenza B/Leningrad/86/93, and 3) Influenza B/Victoria/504/00 stained with ARG10632 anti-Influenza B Haemagglutinin antibody [InB190] at 3 µg/ml dilution.



ARG10632 anti-Influenza B Haemagglutinin antibody [InB190] I-ELISA image

Indirect ELISA: Titration curves of ARG10632 anti-Influenza B Haemagglutinin antibody [InB190] specific to HA of Influenza B virus. Antigens (0.5 µg/well): Influenza B/Tokyo/53/99, Influenza B/Leningrad/86/93, Influenza B/Victoria/504/00, Influenza A: mixture of two strains - A/Shangdong/9/93 and A/New Caledonia/20/99.