

**ARG10932**  
**anti-EBV LMP1 antibody [CS 1-4]**Package: 100 µg  
Store at: -20°C

### Summary

|                     |                                                        |
|---------------------|--------------------------------------------------------|
| Product Description | Mouse Monoclonal antibody [CS 1-4] recognizes EBV LMP1 |
| Tested Reactivity   | EBV                                                    |
| Tested Application  | IHC-P, IP, WB                                          |
| Host                | Mouse                                                  |
| Clonality           | Monoclonal                                             |
| Clone               | CS 1-4                                                 |
| Isotype             | IgG1                                                   |
| Target Name         | EBV LMP1                                               |
| Species             | Virus                                                  |
| Immunogen           | Epstein-Barr Virus, Latent Membrane Protein 1          |
| Conjugation         | Un-conjugated                                          |

### Application Instructions

|                   |                                                                                                                                            |                 |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| Application table | Application                                                                                                                                | Dilution        |
|                   | IHC-P                                                                                                                                      | 1:20            |
|                   | IP                                                                                                                                         | Assay-dependent |
|                   | WB                                                                                                                                         | 1:200           |
| Application Note  | * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist. |                 |

### Properties

|                     |                                                                                                                                                                                                                                                                                                                            |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form                | Liquid                                                                                                                                                                                                                                                                                                                     |
| Purification        | Purified by affinity chromatography.                                                                                                                                                                                                                                                                                       |
| Buffer              | PBS and 0.02% Sodium azide.                                                                                                                                                                                                                                                                                                |
| Preservative        | 0.02% Sodium azide                                                                                                                                                                                                                                                                                                         |
| Storage instruction | For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use. |
| Note                | For laboratory research only, not for drug, diagnostic or other use.                                                                                                                                                                                                                                                       |

ARG10932 anti-EBV LMP1 antibody [CS 1-4] WB image

Western blot: LMP1 expression plasmid transfected 293 cells (left) or vector control 293 cells (right). Cell lysates were stained with ARG10932 anti-EBV LMP1 antibody [CS 1-4].

