

ARG11160 anti-Rabies virus Glycoprotein antibody [4F1]

Package: 200 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [4F1] recognizes Rabies virus Glycoprotein
Tested Reactivity	Virus
Tested Application	EIA, ELISA, ICC/IF, WB
Specificity	MAbs 4F1 are specific to Rabies virus glycoprotein.
Host	Mouse
Clonality	Monoclonal
Clone	4F1
Isotype	IgG2b
Target Name	Rabies virus Glycoprotein
Immunogen	purified rabies virus
Conjugation	Un-conjugated

Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Protein A affinity purified.
Buffer	PBS (pH 7.4) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	1.0-2.0 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Microbiology and Infectious Disease antibody
---------------	--