

ARG20610 Mouse IgG2a Isotype Control antibody [HOPC-1] (APC-Cyanine 5.5)

Package: 50 tests

Store at: 4°C

Summary

Product Description	APC-Cyanine 5.5-conjugated Mouse Monoclonal antibody [HOPC-1] as a negative control antibody for Mouse IgG2a
Tested Application	ELISA, FACS, FLISA, ICC/IF, IHC-Fr, WB
Host	Mouse
Clonality	Monoclonal
Clone	HOPC-1
Isotype	IgG2a, lambda
Target Name	Mouse IgG2a
Conjugation	APC-Cyanine 5.5

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	10 µl/10 ⁶ cells
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.