

ARG20619 Mouse IgG2b Isotype Control antibody [A-1] (PE-Cyanine 7)

Package: 50 tests
Store at: 4°C

Summary

Product Description	PE-Cyanine 7-conjugated Mouse Monoclonal antibody [A-1] as a negative control antibody for Mouse IgG2b
Tested Application	ELISA, FACS, FLISA, IHC-Fr
Specificity	Chicken IgA
Host	Mouse
Clonality	Monoclonal
Clone	A-1
Isotype	IgG2b, kappa
Target Name	Mouse IgG2b
Target Ig	Chicken IgA
Conjugation	PE-Cyanine 7

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	10 µl/10 ⁶ cells
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.