

## ARG20644 Mouse IgG Isotype Control antibody (APC)

Package: 50 µg  
Store at: 4°C

### Summary

Product Description	APC-conjugated Mouse Polyclonal antibody as a negative control antibody for Mouse IgG
Tested Application	BL, FACS, FLISA, IHC-P, WB
Host	Mouse
Clonality	Polyclonal
Isotype	IgG
Target Name	Mouse IgG
Conjugation	APC


### Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	FACS	< 0.1 µg/10 <sup>6</sup> cells
	FLISA	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG20644 Mouse IgG Isotype Control antibody (APC) FLISA image

FLISA: The plate was coated with Goat anti-Mouse IgG antibody, [ARG21539](#) Goat anti-Mouse IgA antibody (pre-adsorbed) and [ARG21517](#) Goat anti-Mouse IgM antibody (pre-adsorbed). Serially diluted ARG20644 Mouse IgG Isotype Control antibody (APC) was captured and fluorescence intensity quantified.