

ARG20703 Rat IgG2a Isotype Control antibody [KLH/G2a-1-1] (Biotin)

Package: 100 µg
Store at: 4°C

Summary

Product Description	Biotin-conjugated Rat Monoclonal antibody [KLH/G2a-1-1] as a negative control antibody for Rat IgG2a
Tested Application	BL, ELISA, FACS, FLISA, IHC-Fr, IHC-P, IP
Specificity	KLH
Host	Rat
Clonality	Monoclonal
Clone	KLH/G2a-1-1
Isotype	IgG2a, kappa
Target Name	Rat IgG2a
Conjugation	Biotin

Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	ELISA	Assay-dependent
	FACS	< 1 µg/10 ⁶ cells
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	IP	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.