

ARG20706
Rat IgG2a Isotype Control antibody [KLH/G2a-1-1] (low endotoxin)Package: 100 µg
Store at: -20°C

Summary

Product Description	Azide free and low endotoxin Rat Monoclonal antibody [KLH/G2a-1-1] as a negative control antibody for Rat IgG2a
Tested Application	BL, ELISA, FACS, FLISA, IHC-Fr, IHC-P, IP
Specificity	KLH
Host	Rat
Clonality	Monoclonal
Clone	KLH/G2a-1-1
Isotype	IgG2a, kappa
Target Name	Rat IgG2a
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	ELISA	Assay-dependent
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	IP	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Purified
Purification Note	Low endotoxin
Buffer	PBS
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.