

ARG20743 Syrian Hamster IgG Isotype Control antibody (APC)

Package: 50 µg
Store at: 4°C

Summary

Product Description	APC-conjugated Syrian Hamster Polyclonal antibody as a negative control antibody for Hamster IgG
Tested Application	ELISA, FACS, FLISA
Host	Syrian Hamster
Clonality	Polyclonal
Isotype	IgG
Target Name	Syrian Hamster IgG
Conjugation	APC

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	< 0.2 µg/10 ⁶ cells
	FLISA	Assay-dependent
	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Application Note		

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.