

ARG20745 Syrian Hamster IgG Isotype Control antibody (Cyanine 5)

Package: 50 µg
Store at: 4°C

Summary

| | |
|---------------------|--|
| Product Description | Cyanine 5-conjugated Syrian Hamster Polyclonal antibody as a negative control antibody for Hamster IgG |
| Tested Application | ELISA, FACS, FLISA |
| Host | Syrian Hamster |
| Clonality | Polyclonal |
| Isotype | IgG |
| Target Name | Syrian Hamster IgG |
| Conjugation | Cyanine 5 |

Application Instructions

| Application table | Application | Dilution |
|-------------------|--|--------------------------------|
| | ELISA | Assay-dependent |
| | FACS | < 0.2 µg/10 ⁶ cells |
| | FLISA | Assay-dependent |
| Application Note | * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist. | |

Properties

| | |
|---------------------|--|
| Form | Liquid |
| Buffer | PBS and 0.1% Sodium azide. |
| Preservative | 0.1% Sodium azide |
| Concentration | 0.1 mg/ml |
| Storage instruction | Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |