

## ARG21185 anti-CD28 antibody [AV7] (PE)

Package: 50 µg  
Store at: 4°C

### Summary

Product Description	PE-conjugated Mouse Monoclonal antibody [AV7] recognizes CD28
Tested Reactivity	Chk, Turkey
Tested Application	FACS, ICC/IF, IHC-Fr
Specificity	Chicken/Turkey CD28. The clone AV7 also reacts with turkey CD28 as demonstrated by flow cytometry.
Host	Mouse
Clonality	Monoclonal
Clone	AV7
Isotype	IgG1, kappa
Target Name	CD28
Species	Chicken
Immunogen	Chicken peripheral blood lymphocytes
Conjugation	PE
Alternate Names	Tp44; TP44; T-cell-specific surface glycoprotein CD28; CD antigen CD28

### Application Instructions

Application table	Application	Dilution
	FACS	< 0.2 µg/10 <sup>6</sup> cells
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.


## Bioinformation

---

Database links	<a href="#">GeneID: 396249 Chicken</a> <a href="#">Swiss-port # P31043 Chicken</a>
Gene Symbol	CD28
Gene Full Name	CD28 molecule
Background	The protein encoded by this gene is essential for T-cell proliferation and survival, cytokine production, and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jul 2011]
Function	Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival. [UniProt]
Research Area	Developmental Biology antibody; Immune System antibody; Cytotoxic T Cell Surface Study antibody
Calculated Mw	25 kDa
PTM	CD40LG induces tyrosine phosphorylation of isoform 3.

## Images

---



ARG21185 anti-CD28 antibody [AV7] (PE) FACS image

Flow Cytometry: Chicken peripheral blood mononuclear cells stained with [ARG21185](#) anti-CD28 antibody [AV7] (PE) and [ARG21151](#) anti-CD3 epsilon antibody [CT-3] (FITC).