

ARG21420 anti-IL10 antibody [JES3-9D7] (low endotoxin)

Package: 100 µg
Store at: -20°C


Summary

Product Description	Azide free and low endotoxin Rat Monoclonal antibody [JES3-9D7] recognizes IL10
Tested Reactivity	Hu, R. Mk, Virus
Tested Application	ELISA, ELISPOT, FACS, IHC-Fr, IHC-P, Neut, WB
Specificity	IL10 antibody clone JES3-9D7 reacts to Human/Viral/Rhesus/Cynomolgus/Sooty Mangabey IL-10, the antibody also reacts to vIL-10 from Epstein-Barr virus (EBV) BCRF1 gene.
Host	Rat
Clonality	Monoclonal
Clone	JES3-9D7
Isotype	IgG1, kappa
Target Name	IL10
Species	Human
Immunogen	COS-expressed human IL-10
Conjugation	Un-conjugated
Alternate Names	IL10A; GVHDS; TGIF; IL-10; Cytokine synthesis inhibitory factor; CSIF; Interleukin-10

Application Instructions

Application table	Application	Dilution
	ELISA	< 5 µg/ml
	ELISPOT	Assay-dependent
	FACS	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	Neut	Assay-dependent
	WB	Assay-dependent

Application Note



ARG21420 anti-IL10 antibody [JES3-9D7] (low endotoxin) standard curve image

ELISA: Human IL-10 detected by [ARG21419](#) anti-IL10 antibody [JES3-9D7] as capture antibody ([ARG21420](#): Low endotoxin version), and [ARG21422](#) anti-IL10 antibody [JES3-12G8] (Biotin) as detection antibody, follow by incubation with streptavidin-HRP.