

## ARG21513 Rat anti-Mouse IgE antibody [23G3] (HRP)

Package: 500 µl  
Store at: 4°C

### Summary

Product Description	HRP-conjugated Rat Monoclonal antibody [23G3] recognizes Mouse IgE
Tested Reactivity	Ms
Tested Application	ELISA, ELISPOT, FACS, FLISA, IHC-Fr, IHC-Wmt, Neut, WB
Specificity	Mouse IgE
Host	Rat
Clonality	Monoclonal
Clone	23G3
Isotype	IgG1, kappa
Target Name	IgE
Species	Mouse
Immunogen	Monoclonal anti-DNP mouse IgE
Conjugation	HRP

### Application Instructions

Application table	Application	Dilution
	ELISA	1:1000 - 1:2000
	ELISPOT	Assay-dependent
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-Wmt	Assay-dependent
	Neut	Assay-dependent
	WB	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	50% PBS (pH 7.4) and 50% Glycerol
Stabilizer	50% Glycerol

**Storage instruction** Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

**Note** For laboratory research only, not for drug, diagnostic or other use.