

ARG21618 Rat anti-Mouse IgA antibody [11-44-2] (AP)

Package: 500 µl
Store at: 4°C

Summary

Product Description	AP-conjugated Rat Monoclonal antibody [11-44-2] recognizes Mouse IgA
Tested Reactivity	Ms
Tested Application	ELISA, FACS, FLISA, IHC-Fr, IHC-P, WB
Host	Rat
Clonality	Monoclonal
Clone	11-44-2
Isotype	IgG1, kappa
Target Name	IgA
Species	Mouse
Conjugation	AP


Application Instructions

Application table	Application	Dilution
	ELISA	1:1000 - 1:2000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl ₂ and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG21618 Rat anti-Mouse IgA antibody [11-44-2] (AP) ELISA image

ELISA: The plate was coated with purified Mouse IgM, IgG, and IgA. Immunoglobulins were detected with serially diluted ARG21618 Rat anti-Mouse IgA antibody [11-44-2] (AP).