

ARG21625 Rat anti-Mouse Kappa Light Chain antibody [187.1] (HRP)

Package: 500 µl
Store at: 4°C

Summary


Product Description	HRP-conjugated Rat Monoclonal antibody [187.1] recognizes Mouse Kappa Light Chain
Tested Reactivity	Ms
Tested Application	ELISA, FACS, FLISA, IHC-Fr, Puri, WB
Host	Rat
Clonality	Monoclonal
Clone	187.1
Isotype	IgG1, kappa
Target Name	Kappa Light Chain
Species	Mouse
Immunogen	MPC-11 immunoglobulin
Conjugation	HRP

Application Instructions

Application table	Application	Dilution
	ELISA	1:4000 - 1:8000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	Puri	Assay-dependent
	WB	Assay-dependent
Application Note	WB: Under non-reducing condition. * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	50% PBS (pH 7.4) and 50% Glycerol
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG21625 Rat anti-Mouse Kappa Light Chain antibody [187.1] (HRP)
ELISA image

ELISA: The plate was coated with purified Mouse IgGκ, IgMκ, IgAκ, IgGλ, IgMλ, and IgAλ. Immunoglobulins were detected with serially diluted ARG21625 Rat anti-Mouse Kappa Light Chain antibody [187.1] (HRP).