

ARG21699 Mouse anti-Rat IgG2a antibody [2A8F4]

Package: 250 µg
Store at: -20°C

Summary

| | |
|---------------------|--|
| Product Description | Mouse Monoclonal antibody [2A8F4] recognizes Rat IgG2a |
| Tested Reactivity | Rat |
| Tested Application | ELISA, ELISPOT, FACS, FLISA, IHC-Fr |
| Host | Mouse |
| Clonality | Monoclonal |
| Clone | 2A8F4 |
| Isotype | IgG1, kappa |
| Target Name | IgG2a |
| Species | Rat |
| Immunogen | Rat IgG2a hybridoma |
| Conjugation | Un-conjugated |

Application Instructions

| Application table | Application | Dilution |
|-------------------|-------------|-----------------|
| | ELISA | Assay-dependent |
| | ELISPOT | Assay-dependent |
| | FACS | Assay-dependent |
| | FLISA | Assay-dependent |
| | IHC-Fr | Assay-dependent |

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

| | |
|---------------------|---|
| Form | Liquid |
| Buffer | BBS (pH 8.2) |
| Concentration | 0.5 mg/ml |
| Storage instruction | For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |