

## ARG21720 Goat anti-Rabbit IgM antibody (AP), pre-adsorbed

Package: 500 µl  
Store at: 4°C

### Summary

Product Description	AP-conjugated Goat Polyclonal antibody recognizes Rabbit IgM
Tested Reactivity	Rb
Tested Application	ELISA, ELISPOT, FACS, FLISA, IHC-Fr, IHC-P, WB
Specificity	The antibody reacts with the heavy chain of Rabbit IgM. The antibody is pre-adsorbed with Rabbit IgG, so the antibody may not react with Rabbit IgG, but may react with immunoglobulins from other species.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgM
Species	Rabbit
Conjugation	AP

### Application Instructions

Pre Adsorbed	Rabbit IgG.	
Application table	Application	Dilution
	ELISA	1:2000 - 1:4000
	ELISPOT	Assay-dependent
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl <sub>2</sub> and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol

Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images

