

ARG21832 Mouse anti-Human IgG2 (Fd) antibody [HP6014]

Package: 250 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [HP6014] recognizes Human IgG2 (Fd)
Tested Reactivity	Hu
Tested Application	Depletion, ELISA, FACS, FLISA, IHC-Fr, IHC-P, IP, Puri, WB
Host	Mouse
Clonality	Monoclonal
Clone	HP6014
Isotype	IgG1, kappa
Target Name	IgG2 (Fd)
Species	Human
Immunogen	Human IgG2 myeloma protein
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	Depletion	Assay-dependent
	ELISA	Assay-dependent
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	IP	Assay-dependent
	Puri	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot

and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.