

ARG21834 Mouse anti-Human IgG2 (Fd) antibody [HP6014] (AP)

Package: 500 µl
Store at: 4°C

Summary

Product Description	AP-conjugated Mouse Monoclonal antibody [HP6014] recognizes Human IgG2 (Fd)
Tested Reactivity	Hu
Tested Application	Depletion, ELISA, FACS, FLISA, IHC-Fr, IHC-P, Puri, WB
Host	Mouse
Clonality	Monoclonal
Clone	HP6014
Isotype	IgG1, kappa
Target Name	IgG2 (Fd)
Species	Human
Immunogen	Human IgG2 myeloma protein
Conjugation	AP

Application Instructions

Application table	Application	Dilution
	Depletion	Assay-dependent
	ELISA	1:500 - 1:1000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	Puri	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl ₂ and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid

repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.