

## ARG21845 Mouse anti-Human IgA2 (Fc) antibody [A9604D2] (FITC)

Package: 250 µg  
Store at: 4°C

### Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [A9604D2] recognizes Human IgA2 (Fc)
Tested Reactivity	Hu
Tested Application	Depletion, FACS, FLISA, IHC-Fr, IHC-P, WB
Host	Mouse
Clonality	Monoclonal
Clone	A9604D2
Isotype	IgG1, kappa
Target Name	IgA2 (Fc)
Species	Human
Immunogen	Purified human IgA2m(1) myeloma protein
Conjugation	FITC

### Application Instructions


Application table	Application	Dilution
	Depletion	Assay-dependent
	FACS	Assay-dependent
	FLISA	1:200 - 1:400
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

## Images



ARG21845 Mouse anti-Human IgA2 (Fc) antibody [A9604D2] (FITC)  
FLISA image

FLISA: The plate was coated with purified Human IgA1, IgA2, IgD, IgG and IgM. Immunoglobulins were detected with serially diluted ARG21845 Mouse anti-Human IgA2 (Fc) antibody [A9604D2] (FITC).