

## ARG21931 Goat anti-Hamster IgG (H+L) antibody (HRP)

Package: 500 µl  
Store at: 4°C

### Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes Hamster IgG (H+L)
Tested Reactivity	Hm
Tested Application	ELISA, FLISA, WB
Specificity	Reacts with the heavy and light chains of hamster IgG; may only react with Syrian hamster IgG
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Hamster
Conjugation	HRP

### Application Instructions

Application table	Application	Dilution
	ELISA	1:4000 - 1:8000
	FLISA	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	50% PBS (pH 7.4) and 50% Glycerol
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.