

## ARG21958 Goat anti-Human Kappa Light Chain antibody (HRP), pre-adsorbed

Package: 500 µl  
Store at: 4°C

### Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes Human Kappa Light Chain
Tested Reactivity	Hu
Tested Application	ELISA, ELISPOT, FACS, FLISA, ICC/IF, IHC-P, Puri, WB
Specificity	Reacts with human κ light chains. May react with κ light chains from other species.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	Kappa Light Chain
Species	Human
Immunogen	Human Kappa light chains
Conjugation	HRP

### Application Instructions

Pre Adsorbed Human λ light chains

#### Application table

Application	Dilution
ELISA	1:4000 - 1:8000
ELISPOT	Assay-dependent
FACS	Assay-dependent
FLISA	Assay-dependent
ICC/IF	Assay-dependent
IHC-P	Assay-dependent
Puri	Assay-dependent
WB	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	50% PBS (pH 7.4) and 50% Glycerol
Stabilizer	50% Glycerol


**Storage instruction**

Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

**Note**

For laboratory research only, not for drug, diagnostic or other use.

**Images**



ARG21958 Goat anti-Human Kappa Light Chain antibody (HRP) (pre-adsorbed) ELISA image

ELISA: The plate was coated with purified Human IgGκ, IgMκ, IgAκ, IgGλ, IgMλ, and IgAλ. Immunoglobulins were detected with serially diluted ARG21958 Goat anti-Human Kappa Light Chain antibody (HRP) (pre-adsorbed).