

ARG22132 anti-FITC antibody (Biotin)

Package: 250 µg
Store at: 4°C

Summary

Product Description	Biotin-conjugated Sheep Polyclonal antibody recognizes FITC
Tested Reactivity	Chem
Tested Application	ELISA, FACS, FLISA, IHC-Fr
Specificity	Reacts with fluorescein isothiocyanate (FITC)
Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Target Name	FITC
Species	Others
Immunogen	Fluorescein isothiocyanate (Isomer I)
Conjugation	Biotin

Application Instructions

Application table	Application	Dilution
	ELISA	1:5000 - 1:20000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.