

## ARG22326 anti-Ly49I antibody [YLI-90]

Package: 100 µg  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [YLI-90] recognizes Ly49I
Tested Reactivity	Ms
Tested Application	FACS
Specificity	Mouse Ly-49I
Host	Mouse
Clonality	Monoclonal
Clone	YLI-90
Isotype	IgG1, kappa
Target Name	Ly49I
Species	Mouse
Immunogen	CHO-K1 cells transfected with B6 allele of Ly49I gene
Conjugation	Un-conjugated
Alternate Names	Ly49I; LY49I1

### Application Instructions

Application table	Application	Dilution
	FACS	< 1 ug/10 <sup>6</sup> cells
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Gene Symbol	Klra9
-------------	-------

Gene Full Name	killer cell lectin-like receptor subfamily A, member 9
Function	Receptor on natural killer (NK) cells for class I MHC. [UniProt]
Calculated Mw	31 kDa