

ARG22402
anti-VZV IE62 antibody [8B11]Package: 50 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [8B11] recognizes VZV IE62
Tested Reactivity	Virus
Tested Application	ICC/IF
Host	Mouse
Clonality	Monoclonal
Clone	8B11
Isotype	IgG1, kappa
Target Name	VZV IE62
Species	Virus
Immunogen	Varicella zoster virus
Conjugation	Un-conjugated


Application Instructions

Application table	Application	Dilution
	ICC/IF	1:200

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.


Properties

Form	Liquid
Purification	Affinity purified.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG22402 anti-VZV IE62 antibody [8B11] ICC/IF image

Immunofluorescence: VZV infected Human fibroblasts stained with ARG22402 anti-VZV IE62 antibody [8B11] (green).



ARG22402 anti-VZV IE62 antibody [8B11] ICC/IF image

Immunofluorescence: HSV-1, HSV-2, CMV, EBV, VZV, and HHV-6 infected cells stained with ARG22402 anti-VZV IE62 antibody [8B11] at 1:100 dilution.
