

ARG22464 anti-CLEC7A antibody [2A11] (Biotin)

Package: 50 µg
Store at: 4°C

Summary

Product Description	Biotin-conjugated Rat Monoclonal antibody [2A11] recognizes CLEC7A This antibody recognizes murine beta-glucan receptor, also known as Dectin-1. Dectin-1 is predominantly expressed by cells of the monocyte/macrophage and neutrophil lineages, but also at lower levels by dendritic cells and a subpopulation of T cells. As a major leucocyte receptor for beta-glucan this molecule may have a key role in the immunomodulatory effects of beta-glucans and in the host response to fungal pathogens. Dectin-1 may stimulate reactive oxygen production in macrophages via the protein tyrosine kinase known as Syk. Rat anti Mouse Dectin-1 antibody, clone 2A11 inhibits the binding of zymosan to macrophages via the beta-glucan receptor.
Tested Reactivity	Ms
Tested Application	FACS
Host	Rat
Clonality	Monoclonal
Clone	2A11
Isotype	IgG2b
Target Name	CLEC7A
Species	Mouse
Immunogen	Dectin-1 transfected NIH3T3 cells and recombinant soluble Dectin-1.
Conjugation	Biotin
Alternate Names	C-type lectin superfamily member 12; CD369; Dectin-1; Dendritic cell-associated C-type lectin 1; CLECSF12; CANDF4; BGR; Beta-glucan receptor; C-type lectin domain family 7 member A; DC-associated C-type lectin 1; DECTIN1

Application Instructions

Application table	<table> <tr> <th>Application</th><th>Dilution</th></tr> <tr> <td>FACS</td><td>1:10 - 1:50</td></tr> </table>	Application	Dilution	FACS	1:10 - 1:50
Application	Dilution				
FACS	1:10 - 1:50				
Application Note	FACS: Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.				

Properties


Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS, 0.09% Sodium azide and 1% BSA.
Preservative	0.09% Sodium azide
Stabilizer	1% BSA

Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	Clec7a
Gene Full Name	C-type lectin domain family 7, member a
Background	This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. The encoded glycoprotein is a small type II membrane receptor with an extracellular C-type lectin-like domain fold and a cytoplasmic domain with an immunoreceptor tyrosine-based activation motif. It functions as a pattern-recognition receptor that recognizes a variety of beta-1,3-linked and beta-1,6-linked glucans from fungi and plants, and in this way plays a role in innate immune response. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. This gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region. [provided by RefSeq, Jul 2008]
Function	Lectin that functions as pattern receptor specific for beta-1,3-linked and beta-1,6-linked glucans, such as cell wall constituents from pathogenic bacteria and fungi. Necessary for the TLR2-mediated inflammatory response and for TLR2-mediated activation of NF-kappa-B. Enhances cytokine production in macrophages and dendritic cells. Mediates production of reactive oxygen species in the cell. Mediates phagocytosis of C.albicans conidia. Binds T-cells in a way that does not involve their surface glycans and plays a role in T-cell activation. Stimulates T-cell proliferation (By similarity). [UniProt]
Calculated Mw	28 kDa
PTM	Phosphorylated on tyrosine residues in response to glucan binding.

Images



ARG22464 anti-CLEC7A antibody [2A11] (Biotin) FACS image

Flow Cytometry: Mouse peripheral blood granulocytes stained with ARG22464 anti-CLEC7A antibody [2A11] (Biotin).