

ARG22923 anti-V5 tag antibody [SV5-Pk1] (Biotin)

Package: 50 µg
Store at: 4°C

Summary

Product Description	Biotin-conjugated Mouse Monoclonal antibody [SV5-Pk1] recognizes V5 tag V5-Tag antibody, clone SV5-Pk1 recognizes the sequence, IPNPLLGL, present on the P/V proteins of the paramyxovirus, SV5. It detects recombinant proteins, including transmembrane and secreted proteins tagged with this sequence.
Tested Reactivity	Other
Tested Application	ELISA, IHC-Fr, WB
Host	Mouse
Clonality	Monoclonal
Clone	SV5-Pk1
Isotype	IgG2a
Target Name	V5 tag
Species	Others
Immunogen	Paramyxovirus Simian-Virus 5 (SV5)
Conjugation	Biotin

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	IHC-Fr	Assay-dependent
	WB	1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS and 0.09% Sodium azide
Preservative	0.09% Sodium azide
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.