

ARG23213 anti-hCG (alpha 5 epitope) antibody [INN-hFSH-158]

Package: 250 µg

Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [INN-hFSH-158] recognizes hCG (alpha 5 epitope) Mouse anti hCG (alpha 5 Epitope) antibody, clone INN-hFSH-158 recognizes the alpha 5 epitope on the alpha subunit of hCG, hLH, hFSH, hISH and the free alpha subunit. The alpha 5 epitope is destroyed upon iodination.
Tested Reactivity	Hu
Tested Application	ELISA, RIA
Host	Mouse
Clonality	Monoclonal
Clone	INN-hFSH-158
Isotype	IgG
Target Name	hCG (alpha 5 epitope)
Species	Human
Immunogen	Human FSH.
Conjugation	Un-conjugated
Alternate Names	FSH-alpha; LSH-alpha; FSHA; Thyrotropin alpha chain; LHA; CG-alpha; GPHA1; TSHA; Chorionic gonadotrophin subunit alpha; Luteinizing hormone alpha chain; TSH-alpha; Choriogonadotropin alpha chain; GPHA; CG-ALPHA; Anterior pituitary glycoprotein hormones common subunit alpha; Follicle-stimulating hormone alpha chain; HCG; Glycoprotein hormones alpha chain; Thyroid-stimulating hormone alpha chain; Follicle-stimulating hormone alpha chain; Lutropin alpha chain

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	RIA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purification with Protein A.
Buffer	PBS and 0.09% Sodium azide.
Preservative	0.09% Sodium azide
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot

and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	CGA
Gene Full Name	glycoprotein hormones, alpha polypeptide
Background	The four human glycoprotein hormones chorionic gonadotropin (CG), luteinizing hormone (LH), follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH) are dimers consisting of alpha and beta subunits that are associated noncovalently. The alpha subunits of these hormones are identical, however, their beta chains are unique and confer biological specificity. The protein encoded by this gene is the alpha subunit and belongs to the glycoprotein hormones alpha chain family. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011]
Calculated Mw	13 kDa
Cellular Localization	Secreted. [UniProt]