

## ARG23438 anti-CD226 antibody [DX11] (low endotoxin)

Package: 250 µg  
Store at: -20°C

### Summary

Product Description	Azide free and low endotoxin Mouse Monoclonal antibody [DX11] recognizes CD226 Mouse anti Human CD226 antibody, clone DX11 recognizes human CD226, a ~65 kDa glycoprotein, also known as DNAM1 (DNAX accessory molecule-1). CD226 is broadly expressed on T-cells, NK cells, platelets, monocytes and a subset of B cells. CD226 is also expressed by a subset of CD3 positive thymocytes. Mouse anti Human CD226 antibody, clone DX11 is reported to inhibit T- and NK cell mediated cytotoxicity against tumor cell targets and to block TNF alpha and IFN gamma secretion by alloantigen-specific T-cells (Kojima et al. 2003).
Tested Reactivity	Hu
Tested Application	FACS, FuncSt, IP
Host	Mouse
Clonality	Monoclonal
Clone	DX11
Isotype	IgG1
Target Name	CD226
Species	Human
Immunogen	Human cytotoxic T lymphocyte clone.
Conjugation	Un-conjugated
Alternate Names	TLISA1; CD antigen CD226; PTA1; CD226 antigen; DNAM-1; DNAX accessory molecule 1; DNAM1

### Application Instructions

Application table	<table><thead><tr><th>Application</th><th>Dilution</th></tr></thead><tbody><tr><td>FACS</td><td>1:50 - 1:200</td></tr><tr><td>FuncSt</td><td>Assay-dependent</td></tr><tr><td>IP</td><td>Assay-dependent</td></tr></tbody></table>	Application	Dilution	FACS	1:50 - 1:200	FuncSt	Assay-dependent	IP	Assay-dependent
Application	Dilution								
FACS	1:50 - 1:200								
FuncSt	Assay-dependent								
IP	Assay-dependent								
Application Note	FACS: Use 10 µl of the suggested working dilution to label 10 <sup>6</sup> cells in 100 µl. * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.								

### Properties

Form	Liquid
Purification	Purification with Protein G.
Purification Note	Low endotoxin
Buffer	PBS
Concentration	1 mg/ml

<b>Storage instruction</b>	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
<b>Note</b>	For laboratory research only, not for drug, diagnostic or other use.

## Bioinformation

---

<b>Gene Symbol</b>	CD226
<b>Gene Full Name</b>	CD226 molecule
<b>Background</b>	This gene encodes a glycoprotein expressed on the surface of NK cells, platelets, monocytes and a subset of T cells. It is a member of the Ig-superfamily containing 2 Ig-like domains of the V-set. The protein mediates cellular adhesion of platelets and megakaryocytic cells to vascular endothelial cells. The protein also plays a role in megakaryocytic cell maturation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2015]
<b>Function</b>	Receptor involved in intercellular adhesion, lymphocyte signaling, cytotoxicity and lymphokine secretion mediated by cytotoxic T-lymphocyte (CTL) and NK cell. [UniProt]
<b>Research Area</b>	Immune System antibody; T Cell Intercellular Adhesion Molecule antibody
<b>Calculated Mw</b>	39 kDa