

ARG23674 anti-Influenza A virus Nucleoprotein N1 antibody [A1]

Package: 50 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [A1] recognizes Influenza A virus Nucleoprotein N1. This product recognizes the Influenza type A nucleoprotein. It demonstrates stronger binding with N1 type influenza A. No cross reactivity is seen with influenza B or other respiratory viruses. Reactivity with virus isolated from other species has not been confirmed with the exception of Avian H5N1 strain.
Tested Reactivity	Influenza A
Tested Application	ICC/IF, IHC-Fr
Host	Mouse
Clonality	Monoclonal
Clone	A1
Isotype	IgG
Target Name	Influenza A virus Nucleoprotein N1
Species	Virus
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified
Buffer	PBS and 0.09% Sodium azide.
Preservative	0.09% Sodium azide
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.