

ARG23694 Goat anti-Horse IgM antibody (HRP)

Package: 250 µg
Store at: -20°C

Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes Horse IgM
Tested Reactivity	Hrs
Tested Application	ELISA, WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgM
Species	Horse
Immunogen	Purified Horse IgM.
Target Ig	IgM
Conjugation	HRP

Application Instructions

Application table	Application	Dilution
	ELISA	1:10000
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified by affinity chromatography.
Buffer	PBS, 0.1% Proclin 300 and 0.2% BSA.
Preservative	0.1% Proclin 300
Stabilizer	0.2% BSA
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.