

ARG23746 anti-H. Pylori antibody

Package: 500 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes H. Pylori.
Tested Reactivity	H. pylori
Tested Application	ELISA, ICC/IF, IHC-P, WB
Specificity	The preparation has not been cross absorbed and may cross react with related microorganisms.
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	H. Pylori
Species	H. pylori
Immunogen	Helicobacter pylori ATCC strain 43504
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	1:5000
	ICC/IF	1:1000
	IHC-P	Assay-dependent
	WB	Assay-dependent
Application Note	WB: Multiple bands between 28 - 87 kDa can be expected. * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	4 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.