

## ARG23754 Rabbit anti-Goat IgG (H+L) antibody (FITC)

Package: 500 µg  
Store at: 4°C

### Summary

|                     |                                                                      |
|---------------------|----------------------------------------------------------------------|
| Product Description | FITC-conjugated Rabbit Polyclonal antibody recognizes Goat IgG (H+L) |
| Tested Reactivity   | Goat                                                                 |
| Tested Application  | ELISA, FLISA, ICC/IF, IHC-Fr, IHC-P, WB                              |
| Specificity         | Reacts with the heavy and light chains of goat IgG                   |
| Host                | Rabbit                                                               |
| Clonality           | Polyclonal                                                           |
| Isotype             | IgG                                                                  |
| Target Name         | IgG (H+L)                                                            |
| Species             | Goat                                                                 |
| Conjugation         | FITC                                                                 |

### Application Instructions

| Application table | Application                                                                                                                                | Dilution        |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
|                   | ELISA                                                                                                                                      | Assay-dependent |
|                   | FLISA                                                                                                                                      | 1:200 - 1:400   |
|                   | ICC/IF                                                                                                                                     | Assay-dependent |
|                   | IHC-Fr                                                                                                                                     | Assay-dependent |
|                   | IHC-P                                                                                                                                      | Assay-dependent |
|                   | WB                                                                                                                                         | Assay-dependent |
|                   | Application Note                                                                                                                           |                 |
|                   | * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist. |                 |

### Properties

|                     |                                                                                                                                                                                                                              |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form                | Liquid                                                                                                                                                                                                                       |
| Purification        | Affinity purification with immunogen.                                                                                                                                                                                        |
| Buffer              | PBS and 0.1% Sodium azide.                                                                                                                                                                                                   |
| Preservative        | 0.1% Sodium azide                                                                                                                                                                                                            |
| Concentration       | 1 mg/ml                                                                                                                                                                                                                      |
| Storage instruction | Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use. |

**Note**

For laboratory research only, not for drug, diagnostic or other use.