

ARG23778 Goat anti-Mouse IgA + IgG + IgM antibody (AP), pre-adsorbed

Package: 500 µl
Store at: 4°C

Summary

Product Description	AP-conjugated Goat Polyclonal antibody recognizes Mouse IgA + IgG + IgM
Tested Reactivity	Ms
Tested Application	ELISA, FACS, FLISA, ICC/IF, IHC-Fr, IHC-P, WB
Specificity	Reacts with the heavy and light chains of mouse IgM, IgG1, IgG2a, IgG2b, IgG3 and IgA as demonstrated by ELISA, FLISA, and/or flow cytometry; minimal cross reactivity with human immunoglobulins
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgA + IgG + IgM
Species	Mouse
Conjugation	AP

Application Instructions


Pre Adsorbed	Human immunoglobulins and pooled sera.	
Application table	Application	Dilution
	ELISA	1:2000 - 1:4000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl ₂ and 50% Glycerol.
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol

Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images



ARG23778 Goat anti-Mouse IgA + IgG + IgM antibody (AP) (pre-adsorbed) IHC-P image

Immunohistochemistry: Paraffin-embedded chicken nasal cavity section infected with H7N1 highly pathogenic avian influenza virus was stained with anti-IAV NP followed by ARG23778 Goat anti-Mouse IgA + IgG + IgM antibody (AP) (pre-adsorbed) and BCIP/NBT and then anti-PGP9.5 followed a secondary antibody and DAB.