

## ARG23831 Mouse anti-Rat Kappa Light Chain antibody [K4F5] (AP)

Package: 500 µl  
Store at: 4°C

### Summary

Product Description	AP-conjugated Mouse Monoclonal antibody [K4F5] recognizes Rat Kappa Light Chain
Tested Reactivity	Rat
Tested Application	ELISA, FACS, FLISA, IHC-Fr, WB
Host	Mouse
Clonality	Monoclonal
Clone	K4F5
Isotype	IgG2a, kappa
Target Name	Kappa Light Chain
Species	Rat
Immunogen	Rat Kappa hybridoma
Conjugation	AP

### Application Instructions

Application table	Application	Dilution
	ELISA	1:1000 - 1:2000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	WB	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Buffer	Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl <sub>2</sub> and 50% Glycerol.
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.