

ARG23886 Goat anti-Human IgA + IgG + IgM antibody

Package: 1 mg
Store at: -20°C

Summary

Product Description	Goat Polyclonal antibody recognizes Human IgA + IgG + IgM
Tested Reactivity	Hu
Tested Application	ELISA, FACS, FLISA, ICC/IF, WB
Specificity	Reacts with the heavy and light chains of Human IgG, IgM, and IgA.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgA + IgG + IgM
Species	Human
Immunogen	Human IgG, IgM, and IgA.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	Assay-dependent
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	BBS (pH 8.2)
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.