

## ARG23905 Goat anti-Human IgG (H+L) antibody, pre-adsorbed

Package: 500 µg  
Store at: -20°C

### Summary

Product Description	Goat Polyclonal antibody recognizes Human IgG (H+L)
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	The antibody reacts with the heavy and light chains of human IgG, and the light chains of human IgM and IgA. The antibody is pre-adsorbed with rabbit, mouse, rat, bovine, horse, hamster, goat, sheep, chicken, and guinea pig serum proteins. The antibody may react with immunoglobulins from other species.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Human
Immunogen	Human IgG
Conjugation	Un-conjugated

### Application Instructions

Pre Adsorbed Rabbit, mouse, rat, bovine, horse, hamster, goat, sheep, chicken, and guinea pig serum proteins.

#### Application table

Application	Dilution
ELISA	< 1 µg/ml

Application Note \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	BBS (pH 8.2)
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.