

ARG23920 Mouse anti-Monkey IgG (H+L) antibody [SB108a]

Package: 250 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [SB108a] recognizes Monkey IgG (H+L)
Tested Reactivity	Mk, R. Mk
Species Does Not React With	Hu, Ms, Rat, Bov, Cat, Chk, Dnk, Goat, Gpig, Hm, Hrs, Pig, Rb, Sheep
Tested Application	ELISA, FACS, FLISA, WB
Specificity	This antibody reacts with Rhesus and cynomolgus IgG. Minimal reactivity to human and rabbit immunoglobulins, goat IgG, and mouse, rat, hamster, guinea pig, sheep, donkey, bovine, horse, porcine, feline, and chicken serum.
Host	Mouse
Clonality	Monoclonal
Clone	SB108a
Isotype	IgG1
Target Name	IgG (H+L)
Immunogen	Rhesus monkey IgG.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	< 1 µg/ml
	FACS	Assay-dependent
	FLISA	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.