

ARG23934 Mouse anti-Human IgG (Fc) antibody [H2] (FITC)

Package: 250 µg
Store at: 4°C

Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [H2] recognizes Human IgG (Fc)
Tested Reactivity	Hu, AGMK, Chimp, R. Mk
Tested Application	ELISA, ELISPOT, FACS, FLISA, WB
Host	Mouse
Clonality	Monoclonal
Clone	H2
Isotype	IgG2b, kappa
Target Name	IgG (Fc)
Species	Human
Immunogen	Human IgG
Conjugation	FITC


Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	ELISPOT	Assay-dependent
	FACS	< 1 µg/10 ⁶ cells
	FLISA	1:200 - 1:600
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG23934 Mouse anti-Human IgG (Fc) antibody [H2] (FITC) FLISA image

FLISA: The plate was coated with purified Human IgG, IgM, and IgA. Immunoglobulins were detected with serially diluted ARG23934 Mouse anti-Human IgG (Fc) antibody [H2] (FITC).