

## ARG24047 anti-Aflatoxin (total) antibody [BCAB1314]

Package: 50 µg  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [BCAB1314] recognizes Aflatoxin (total)
Tested Reactivity	Other
Tested Application	ELISA
Specificity	Cross Reactivity: Aflatoxin B1: 100% Aflatoxin B2: 86% Aflatoxin G1: 101% Aflatoxin G2: 82%
Host	Mouse
Clonality	Monoclonal
Clone	BCAB1314
Isotype	IgG1, lambda
Target Name	Aflatoxin (total)
Species	Others
Immunogen	Aflatoxin (total).
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Purification with Protein G.
Purity	> 95% (by electrophoresis)
Buffer	Neutralized citrate buffer (pH 7.4) and 0.05% Sodium azide.
Preservative	0.05% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.