

ARG24050 anti-Biotin antibody [SB58c]

Package: 100 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [SB58c] recognizes Biotin
Tested Reactivity	Other
Tested Application	ELISA, FACS, FLISA, WB
Host	Mouse
Clonality	Monoclonal
Clone	SB58c
Isotype	IgG1
Target Name	Biotin
Species	Others
Immunogen	KLH-conjugated Biotin
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	1 µg/ml
	FACS	Assay-dependent
	FLISA	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.