

ARG24057 Goat anti-Human IgG (Fd) antibody (FITC), pre-adsorbed

Package: 500 µg
Store at: 4°C

Summary

Product Description	FITC-conjugated Goat Polyclonal antibody recognizes Human IgG (Fd)
Tested Reactivity	Hu
Tested Application	FLISA
Specificity	The antibody reacts with the Fd region of human IgG. The antibody is pre-adsorbed with Human IgG Fc, IgM and IgA, so the antibody may not react with Human IgG Fc, IgM and IgA, but may react with IgG from other species.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (Fd)
Species	Human
Immunogen	Human IgG
Conjugation	FITC

Application Instructions

Pre Adsorbed	Human IgG Fc, IgM and IgA.				
Application table	<table> <tr> <th>Application</th><th>Dilution</th></tr> <tr> <td>FLISA</td><td>1:200 - 1:400</td></tr> </table>	Application	Dilution	FLISA	1:200 - 1:400
Application	Dilution				
FLISA	1:200 - 1:400				
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.				


Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Highlight

Related news:
[Comprehensive anti-Human secondary antibodies](#)

Images



ARG24057 Goat anti-Human IgG (Fd) antibody (FITC) (pre-adsorbed)
FLISA image

FLISA: The plate was coated with purified Human IgG Fab and IgG Fc. Immunoglobulins were detected with serially diluted ARG24057 Goat anti-Human IgG (Fd) antibody (FITC) (pre-adsorbed).